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

Laboratory Equipment of Platinum and Platinum Alloys

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Crucibles and Dishes

In view of the chemical stability and resistance of the precious metals, crucibles and dishes of platinum and platinum alloys are predestined for use as equipment in physical-chemistry investigations.

The most commonly used crucibles, dishes and covers are standardized and are listed in the tables.

To avoid mechanical deformation we deliver the crucibles and dishes in different thicknesses and with a reinforced rim. All dishes can be produced with or without a pouring lip.

Special designs that are required for a wide variety of analytical methods, including standard techniques, also belong to the delivery programme.

We offer special shapes and sizes according to customer's requirements. The decision as to whether platinum or a platinum alloy is better suited to your particular application depends on the aggressiveness of the substances used and on the physical parameters such as temperature, atmosphere, etc. The appendix „Handling Platinum Equipment“ gives the selection criteria regarding the choice of suitable materials.

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Hemispherical dishes having the same capacities as above are also available. Prices on request.

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2) äußere Elektrode / outer electrode

XFA Programme

The demands on elemental analysis have become more severe with the increasing importance of quality assurance. X-ray fluorescence analysis (XFA) has gained a key position because of its high speed and reproducibility.

In order to maintain the required analytical precision, no impurities may be introduced during the preparation of the melted tablets for XFA. This extreme requirement can only be met by the high corrosion resistance of platinum. Added to this is the fact that the molten glass can be poured out of the melting crucible without leaving any deposit and the solidified button can be easily removed from the casting dish. This is achieved by the addition of 5 % gold to the platinum, which reduces wetting. Furthermore the gold addition to the alloy reduces the tendency to recrystallization so that there is no influence of the grain structure on the surface of the melted button.

Electrodes

Platinum-Iridium 90/10 is used as the electrode material to ensure mechanical stability.

Common electrode types

This electrode type consists of two concentric gauze cylinders. The electrode is sand-blasted to enlarge the wire surface area.

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This design consists of a slit gauze cylinder containing a concentrically positioned wire spiral.

Bestell-Nr. Cat. No.	Gesamthöhe Total height [mm]	Elektrodenhöhe Electrode height [mm]	Durchmesser Diameter [mm]	Drahtstärke Wire thickness [mm]	Gewicht ca. approx. weight [g]
EL 01/1 ³⁾	150	50	35	0.12	16.0
EL 01/2 ³⁾	150	50	35	0.25	22.5
EL 02 ⁴⁾	150	70	10	1.1	7.5

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In this electrode the gauze is fixed stably to the central rod by means of a solid cross at top and bottom.

Bestell-Nr. Cat. No.	Gesamthöhe Total height [mm]	Elektrodenhöhe Electrode height [mm]	Durchmesser Diameter [mm]	Drahtstärke Wire thickness [mm]	Gewicht ca. approx. weight [g]
EL 05/1	140	40	20	0.12	11.5
EL 05/2	140	40	20	0.25	14.5
EL 05/1	140	40	30	0.12	20.0
EL 05/2	140	40	30	0.25	26.0
EL 05/1	130	50	45	0.12	35.0
EL 05/2	130	50	45	0.25	45.0

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This type consists of a platinum gauze bent into a U-shape. It is used mainly for the rapid determination of halogens and sulfur.

Bestell-Nr. Cat. No.	Gesamthöhe Total height [mm]	Elektrodenhöhe Electrode height [mm]	Durchmesser Diameter [mm]	Drahtstärke Wire thickness [mm]	Gewicht ca. approx. weight [g]
EL 06	45	35	10	0.76	3

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Other Laboratory Equipment

The main items in this group of products are:

- Calorimetry dishes in different shapes and sizes, including standard DIN designs
- Boats in different sizes with and without handles
- Stainless steel crucible tongs with platinum platelets or platinum shoes at the tips
- Stainless steel forceps with platinum platelets riveted on
- Spatulas in various shapes and lengths.

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Gauzes

Listed below are gauzes of platinum, platinum-rhodium, platinum-iridium and silver which can be used, for example, for catalytic or electrochemical purposes. They are representative of some of the current alloys, mesh sizes and wire thicknesses. Within the usual range of wire diameters, our production programme includes mesh counts ranging from 10 to 3600 meshes per cm².

Legierung Alloy	Maschenzahl Mesh count	Drahtstärke Wire diameter [mm]	Maschenweite Mesh size [mm]	Gewicht ca. approx. weight [g / 100 cm ²]
Pt	16	0.76	1.7	86.0
Pt	1024	0.06	0.25	4.3
Pt	1024	0.076	0.2	6.5
Pt	3600	0.04	0.12	3.6
PtRh 90/10	1024	0.06	0.25	4.0
PtIr 90/10	100	0.25	0.75	22.0
PtIr 90/10	250	0.12	0.5	8.0
Ag	64	0.35	0.9	18.0
Ag	1024	0.12	0.19	8.3
AuPt 10	1024	0.06	0.25	3.9

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Instructions about handling of Platinum equipment are available on request.